

Account: \_\_\_\_\_ Date: \_\_\_\_\_

Farm Name: \_\_\_\_\_ Contact Person: \_\_\_\_\_

Street: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

Phone: \_\_\_\_\_ E-mail: \_\_\_\_\_

**Reorder:** 10 vials x \_\_\_\_\_

### Direction for using the milk sampling kit

1. Order kits by emailing [supply@dairyone.com](mailto:supply@dairyone.com) or check box for re-order. Kits can also be ordered through the online store at [www.dairyone.com](http://www.dairyone.com)
2. Obtain a milk sample for each animal to be tested in one of the pre-served sample vials. To ensure a usable sample each vial should be at least 1/3 full.
3. Close the lid tightly and mark the sample vial with the vial number on top with a permanent marker. Record the ID of the animal with the vial number on the submission form. Record DIM information.
4. Place vials in the provided whirlpak bag. Place whirlpak bag in the bubble envelope and refrigerate samples until ready to ship.
5. Place submission form in the bubble envelope (not in the whirlpak bag) and place label on the envelope.
6. Free shipping for 10 or more samples.

**Send to: Dairy One Animal Health Diagnostic Laboratory  
720 Warren Road, Ithaca, New York 14850**

If you have questions or would like more information please email:  
[dhiarecordservices@dairyone.com](mailto:dhiarecordservices@dairyone.com) | Call: 607-252-7575 | Txt: 607-257-1272

By using this form, the customer confirms that they understand the following:

- A fee of \$15 may be applied if sample submissions require excessive sample handling and/or directions are not followed.
- Days In Milk information is necessary to determine the suitability of samples for testing and to provide proper interpretation of results. By not providing DIM information that conforms to validity criteria, the validity of the test may be affected. Cows  $\geq 5$  DIM.
- If a Sampled Date is not provided, the date of the test will be used.
- Simplified test results in the form of a .csv file may be sent to an unsecure digital area.
- The AHDL uses Mycoplasma bovis DNA PCR Testing method on individual cow samples.
- **Strip samples should be taken aseptically.**
- AHDL uses PCR Testing on individual cow samples to detect the presence of Mycoplasma bovis DNA, a pathogenic bacterial species that can cause bovine mastitis.

Vial #	Animal ID Number	Days in Milk
1		
2		
3		
4		
5		
6		
7		
8		
9		
10		
11		
12		
13		
14		
15		

Vial #	Animal ID Number	Days in Milk
16		
17		
18		
19		
20		
21		
22		
23		
24		
25		
26		
27		
28		
29		
30		

Vial #	Animal ID Number	Days in Milk
31		
32		
33		
34		
35		
36		
37		
38		
39		
40		
41		
42		
43		
44		
45		
46		
47		
48		
49		
50		
51		
52		
53		
54		
55		
56		
57		
58		
59		
60		
61		
62		
63		
64		
65		

Vial #	Animal ID Number	Days in Milk
66		
67		
68		
69		
70		
71		
72		
73		
74		
75		
76		
77		
78		
79		
80		
81		
82		
83		
84		
85		
86		
87		
88		
89		
90		
91		
92		
93		
94		
95		
96		
97		
98		
99		
100		

Herdcode: \_\_\_\_\_